

## Trypan Blue Tissue Culture

---

|                   |                    |                    |                            |
|-------------------|--------------------|--------------------|----------------------------|
| Sample Name:      | <b>Sample Name</b> | Created Date/Time: | <b>2025-04-04 09:00:51</b> |
| Account:          | <b>General</b>     | Passage:           |                            |
| Software Version: | <b>v4.1.8</b>      | Firmware Version:  | <b>11.66</b>               |
| Instrument SN:    | <b>s-07097</b>     | Record:            | <b>20250404090051</b>      |
| Description:      |                    |                    |                            |

## Results

---

|                           |                            |                           |                 |
|---------------------------|----------------------------|---------------------------|-----------------|
| Last Processed Date/Time: | <b>2025-04-04 09:00:51</b> | Gating Min Diameter (um): | <b>6</b>        |
| Gating Max Diameter (um): | <b>17</b>                  | Live Cell Count:          | <b>160</b>      |
| Dead Cell Count:          | <b>275</b>                 | Total Cell Count:         | <b>435</b>      |
| % Viability:              | <b>36.7816</b>             | Mean Diameter (um):       | <b>13.03</b>    |
| Live Cells/mL:            | <b>9.41e+05</b>            | Dead Cells/mL:            | <b>1.62e+06</b> |
| Total Cells/mL:           | <b>2.56e+06</b>            |                           |                 |

## Dilution

---


|                     |            |                         |            |
|---------------------|------------|-------------------------|------------|
| Target Cells/mL:    | <b>N/A</b> | Target Number of Cells: | <b>N/A</b> |
| Target Volume (mL): | <b>0</b>   | Original Volume (mL):   | <b>0</b>   |
| Dilution Guidance:  |            |                         |            |

## Protocol


---

|                             |                         |                             |                |
|-----------------------------|-------------------------|-----------------------------|----------------|
| Protocol:                   | <b>Default Protocol</b> | DirectPipette (or Slide):   | <b>Pipette</b> |
| Chamber Height:             | <b>100</b>              | Dilution Factor:            | <b>2</b>       |
| Protocol Min Diameter (um): | <b>6</b>                | Protocol Max Diameter (um): | <b>17</b>      |
| Live Roundness:             | <b>50</b>               | Dead Roundness:             | <b>25</b>      |
| Stained Threshold:          | <b>35</b>               |                             |                |
| Small Cell Mode:            | <b>False</b>            |                             |                |
| Irregular Cell Mode:        | <b>False</b>            |                             |                |

## Sample Images




Brightfield




Counted


## Sample Histograms




Live Cells (by Size)




Dead Cells (by Size)



Live+Dead Cells (by Size)



Total Cells (by Size)



Cluster Size